

Goat Anti-Mouse IgG H&L(PE)



Cat No.: SJ-SA6680

Product Information

Host

Goat

Reactive Species

Mouse

Applications

IHC-F,IHC-P,ICC/IF,FC,ELISA

Conjugation

PE. Ex:565nm.Em:576nm.

Detailed Information

Recommended Dilution

IHC:1:500-1:1000

IF:1:100-1:1000

FC:1:1000-1:4000

ELISA:1:5000-1:20000

Storage Buffer

Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide

Storage Conditions

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Concentration

N/A

Alternative Names

N/A

NOTICE: For research use only. Not for therapeutic or diagnostic purposes.

Contact

For technical support and original validation data for this product please contact:

T: 400-018-6916

E: support@sanjingbio.com

W: www.sanjingbio.com