

# Rabbit Anti-Mouse IgG H&L (Biotin)



Cat No.: SJ-SA6752

## Product Information

---

Host

Rabbit

Reactive Species

Mouse

Applications

WB,IHC-F,IHC-P,ICC/IF,ELISA

Conjugation

Biotin

## Detailed Information

---

Recommended Dilution

WB:1:5000-1:20000

IF:1:100-1:1000

ELISA:1:5000-1:20000

Storage Buffer

Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.05% Sodium Azide

Storage Conditions

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Concentration

N/A

Alternative Names

N/A

**NOTICE:** For research use only. Not for therapeutic or diagnostic purposes.

## Contact

---

For technical support and original validation data for this product please contact:

T: 400-018-6916

E: support@sanjingbio.com

W: www.sanjingbio.com