

DyLight 550 Conjugated Avidin



Cat No.: SJ-SA6831

Product Information

Host N/A	Reactive Species N/A
Applications IF, ICC/IF, FCM	Conjugation Dylight 550

Detailed Information

Recommended Dilution

IF
ICC/IF
FCM 1:200—800; Optimal working dilutions must be determined by end users.

Storage Buffer

1 mg/ml of DyLight 550 conjugated Avidin; 0.01M PBS, 5 mg/ml BSA, 50% glycerol.

Storage Conditions

At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protected from light.

Concentration

1 mg/ml

Alternative Names

N/A

NOTICE: For research use only. Not for therapeutic or diagnostic purposes.

Contact

For technical support and original validation data for this product please contact:

T: 400-018-6916

E: support@sanjingbio.com

W: www.sanjingbio.com